EN

**Remington  
MEDICAL**

**AUTOMATIC CUTTING NEEDLE**

NAC Needle Family

(NAC, NAC B, NAC M, NAC MLL, NAC U, NAC ULTRA)

SOFT TISSUE BIOPSY NEEDLE

**Indications for Use:**

Automatic Cutting Needles (NAC) are intended to be used by professional medical personal for use in obtaining biopsies from soft tissues such as liver, kidney, prostate, spleen, lymph nodes, and various soft tissue tumors to aid in diagnosis of cancer and/or other tissue abnormalities when used in conjunction with the following biopsy instruments as represented in the table below.

<i>Biopsy Instrument:</i>	<i>Fits with Needle:</i>
Angiotech® Pro-Mag™ 2.2	NAC, NAC B, NAC U
Angiotech® Pro-Mag™ ULTRA	NAC ULTRA
Bard® Biopsy®	NAC, NAC B, NAC U
Bard® MAGNUM®	NAC M, NAC MLL
MBD™-23	NAC, NAC B

**Contraindications:**

1. Intended for all patient populations, however, physician judgment is required when considering a biopsy on patients with bleeding disorders and/or patients receiving anticoagulant medications.
2. Needle should be used by a physician familiar with the possible side effects, typical findings, limitations, indications and contraindications of a core needle biopsy procedure.
3. For single patient use only. Can be utilized at the physician's discretion at the same tissue site to collect multiple tissue samples. Not intended to be reused at multiple tissue sites. If re-used, infection may occur and/or needle set may not correctly fit with applicable biopsy instrument or may bend/break during use.
4. Not intended for use in bone.

**Instructions for Use:**

1. Open the sterile package using proper aseptic technique.
2. Remove the needle from the sterile package and remove the tip protector.
3. Inspect the needle (stylet and cannula) for damage that would prevent proper operation of the needle assembly.
4. Open the lid of the biopsy device and load the needle set into the device.
5. Operate the biopsy device following the manufacturer's operating instructions.
6. If multiple samples are taken, inspect the needle (stylet and cannula) for damage that would prevent proper operation of the needle assembly between each respective tissue sample.

**Warning:** Those who use or may be exposed to needles are at increased risk of needle related injury. Medical devices such as this may present an increased risk of infection and/or needlestick injuries. To mitigate these risks, application of the device by a physician familiar with core needle biopsy is suggested. After use, this product may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable laws and regulations. Handle in a manner which will prevent accidental needle puncture.

Any serious incident which has occurred in relation to the NAC Family should be reported to Remington Medical, Inc. and the competent authority of the Member State in which the user and/or patient is established.

Sterilized by  
Ethylene Oxide

2797

Do Not  
Reuse

Do Not  
Resterilize

EMERGO EUROPE  
Prinsessegracht 20  
2514 AP The Hague  
The Netherlands

EC REP

Do Not Use  
if Package  
is Open or  
Damaged

Rx Only

Manufacturer

Non-Pyrogenic

Remington Medical, Inc.  
6830 Meadowridge Court  
Alpharetta, GA 30005 USA  
+1.800.989.0057  
quality@remmed.com

Angiotech® and Pro-MagM are registered trademarks of Angiotech or its affiliates. Bard®, Biopsy® and MAGNUM® are registered trademarks of C.R. Bard, nc or its affiliates. MBDM is a registered trademark of Promex Technologies, LLC or its affiliates.

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Endast Rx

Remington Medical, Inc.  
6830 Meadowridge Court  
Alpharetta, GA 30005 USA  
+1 800 989.0057  
quality@remmed.com

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Samo na recept

Remington Medical, Inc.  
6830 Meadowridge Court  
Alpharetta, GA 30005 ZDA  
+1.800.989.0057  
quality@remmed.com

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