
Remington[®] MEDICAL

EN

S-101-97 / S-101-97-B / S-201-97 DISPOSABLE EXTENSION CABLE

Basic UDI-DI: 081307902CABLEW7

INTENDED PURPOSE:

The Disposable Extension Cables are defined as surgically non-invasive devices intended for transient use. The Disposable Extension Cable Families are intended to be used by professional medical personnel during pacemaker implantation and temporary pacing. The cables are single-patient use devices and are not intended to be used for more than one patient nor for several procedures as these may increase risk of infection. The S-101-97 / S-101-97-B / S-201-97 Disposable Extension Cable is intended to act as a conduit for the electrical signal generated from a Pacing System Analyzer (PSA) or External Pacing Generator (EPG) to a permanent pacemaker lead or epicardial wire. These cables are used in conjunction with a PSA/EPG unit, permanent pacemaker catheters, epicardial wires or surgical cable adaptors. After connections are completed, temporary pacing can be initiated if connected to an EPG, or if utilizing a PSA, pacing, and sensing thresholds can be measured.

INDICATIONS FOR USE:

The S-101-97 / S-101-97-B / S-201-97 Disposable Extension Cables are electrical extension cables designed to transmit signal from patient-connected electrodes to a PSA, or power or excitation signal from an EPG to patient-connected electrodes. After connections are completed, temporary pacing can be initiated if connected to an EPG, or if utilizing a PSA, pacing and sensing thresholds can be measured. The S-101-97 / S-101-97-B / S-201-97 Disposable Extension Cables do not have any therapeutic or diagnostic indications or claims.

WARNINGS/PRECAUTIONS:

1. Transmittance of micro-shocks from static electricity can be generated through the pacing system. It is recommended that gloves be used when

- handling any part of the pacing system while in use.
2. Damage to cable may occur from improper handling.
 3. DO NOT reuse. This may cause the cables to malfunction and/or to no longer fit with the EPG/PSA unit.
 4. DO NOT use if sterile packaging is damaged or unintentionally opened before use.
 5. DO NOT use for patients weighing less than 6.49 kilograms.

CONTRAINDICATIONS:

The Disposable Extension Cable Family is intended for all patient populations; however, physician judgment is required regarding certain patient populations for pacing procedures (i.e., use in patients with conditions that may affect temporary or permanent pacing).

INSTRUCTIONS FOR USE:

1. With Temporary Pacemaker:

Remove the S-101-97 / S-101-97-B / S-201-97 Disposable Extension Cable from the sterile package and insert the Safe Connect connector into the EPG until it clicks. Then attach the alligator clips to the appropriate temporary pacing wires. NOTE: In dual chamber pacing a second cable is required and connections are made by appropriately connecting the positive (red) and negative (black) leads of each chamber.

2. With Pacing System Analyzer:

Remove the S-101-97 / S-101-97-B / S-201-97 Disposable Extension Cable from the sterile package and insert the Safe-Connect connector into the PSA until it clicks. Then attach the alligator clips to the ends of the permanent pacemaker leads.

a. Single-Chamber Unipolar Leads

Requires one S-101-97 / S-101-97-B / S-201-97 Disposable Extension Cable. Attach the black (negative) clip to the connector pin on the end of lead. Attach the red (positive) clip to the temporary indifferent electrode.

b. Single-Chamber Bipolar Leads

Requires one S-101-97 / S-101-97-B / S-201-97 Disposable Extension Cable. Attach the black (negative) clip to the negative pin on the end of lead. Attach the red (positive) clip to the positive pin or positive ring on the end of lead.

c. Dual-Chamber Unipolar Leads

Requires two S-101-97 / S-101-97-B / S-201-97 Disposable Extension Cables. Attach the black (negative) clip of the atrial cable to the connector pin on the end of the atrial lead and attach the black (negative) clip of the second (ventricular) cable to the connector pin on the end of the ventricular lead. Attach the red (positive) clip of both cables to the temporary indifferent electrode.


d. Dual-Chamber Bipolar Leads

Requires two S-101-97 / S-101-97-B / S-201-97 Disposable Extension Cables. Attach the black (negative) clip of the atrial cable to the connector pin on the end of the atrial lead and attach the black (negative) clip of the second cable to the connector pin on the end of the ventricular lead. Attach the red (positive) clip of atrial cable to the positive pin/ring on the end of the atrial lead and attach the second red (positive) clip of the ventricular cable to the positive pin/ring on the end of the ventricular lead.

DISPOSAL:

1. Dispose of the S-101-97 / S-101-97-B / S-201-97 in an environmentally friendly way as contaminated medical waste.

Any serious incident which has occurred in relation to the S-101-97 / S-101-97-B / S-201-97 Disposable Extension Cable should be reported to Remington Medical, Inc. and to the competent authority of the Member State in which the user and/or patient is established. The Summary of Safety and Clinical Performance (SSCP) for Remington Medical Disposable Extension Cables is available in the European database on medical devices (Eudamed) where it is linked to the Basic UDI-DI. The URL to the Eudamed public website is: <https://ec.europa.eu/tools/eudamed>. The SSCP is also available online at remmed.com

Rx Only

RoHS Compliant
Lead-Free

Remington Medical, Inc.
6830 Meadowridge Court
Alpharetta, GA 30005 USA
+1.800.989.0057
quality@remmed.com

Do Not Re-sterilize

Do Not Re-use

Do Not Use if Packaging is Open or Damaged

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