
Remington[®] MEDICAL

EN

BS-101-97 / BS-101-97-12 DISPOSABLE EXTENSION CABLE

Basic UDI-DI: 081307902CABLEW7

INTENDED PURPOSE:

The Disposable Extension Cables are defined as surgically non-invasive devices intended for transient use. The Disposable Extension Cable Families are intended to be used by professional medical personnel during specific cardiac procedures such as pacemaker implantation and cardiac functionality analysis, etc. The cables are a single-patient use device and are not intended to be used for more than one patient nor for several procedures as these may increase risk of infection. The BS-101-97/BS-101-97-12 Bi-ventricular Disposable Extension Cable is intended to act as a conduit for the electrical signal generated from a Pacing System Analyzer (PSA) to a permanent pacemaker lead. After connections are completed, analyzing, pacing, and sensing thresholds can be measured.

INDICATIONS FOR USE:

The BS-101-97 / BS-101-97-12 Disposable Extension Cables are electrical extension cables designed to transmit signal from patient-connected electrodes to a PSA, or power or excitation signal from an EPG to patient-connected electrodes. After connections are completed, temporary pacing can be initiated if connected to an EPG, or if utilizing a PSA, pacing and sensing thresholds can be measured. The BS-101-97 / BS-101-97-12 Disposable Extension Cable do not have any therapeutic or diagnostic indications or claims.

WARNINGS/PRECAUTIONS:

1. Transmittance of micro-shocks from static electricity can be generated through the pacing system. It is recommended that gloves be used when handling any part of the pacing system while in use.
2. Damage to cable may occur from improper handling.
3. DO NOT reuse. This may cause the cables to malfunction and/or to no longer fit with the PSA unit.
4. DO NOT use if sterile packaging is damaged or unintentionally opened before use.
5. DO NOT use for patients weighing less than 6.49 kilograms.

CONTRAINDICATIONS:

The Disposable Extension Cable Family is intended for all patient populations; however, physician judgment is required regarding certain patient populations for pacing procedures (i.e., use in patients with conditions that may affect temporary or permanent pacing).

INSTRUCTIONS FOR USE:

Remove the BS-101-97 / BS-101-97-12 Disposable Extension Cable from the sterile package and insert the Safe-Connect connector into the PSA or adaptor as follows:

1. Single-Chamber Pacing:

Bi-Ventricular (one unipolar, one bipolar) leads require one BS-101-97 / BS-101-97-12 Disposable Extension Cable. Connect the cable for the ventricular lead system to the ventricular connector block. Attach the alligator clips at the proximal end of the cable in either a unipolar or bipolar configuration.

a. *Single-Chamber Bi-Ventricular (one unipolar, one bipolar)*

Attach one black (negative) clip to the negative pin on the end of the bipolar lead. Attach the second black (negative) clip to the negative connector pin on the end of the unipolar lead. Attach the red (positive) clip to the positive pin or positive ring on the end of the bipolar lead.

2. Dual-Chamber Pacing:

Bi-Ventricular (one unipolar, one bipolar) leads require one BS-101-97 / BS-101-97-12 Disposable Extension Cable for Bi-Ventricular leads, and one (S-101-97 or S-201-97) Disposable Extension Cable for the atrial lead.

a. *Dual-Chamber Bi-Ventricular (one unipolar, one bipolar)*

Attach one black (negative) clip to the negative pin on the end of the bipolar ventricular lead. Attach the second black (negative) clip to the connector pin on the end of the unipolar ventricular lead. Attach the red (positive) clip to the positive pin or positive ring on the end of the atrial lead.

DISPOSAL:

1. Dispose of the BS-101-97 / BS-101-97-12 Disposable Extension Cable in an environmentally friendly way as contaminated medical waste.

Any serious incident which has occurred in relation to the BS-101-97 / BS-101-97-12 should be reported to Remington Medical, Inc. and to the competent authority of the Member State in which the user and/or patient is established. The Summary of Safety and Clinical Performance (SSCP) for Remington Medical Disposable Extension Cables is available in the European database on medical devices (Eudamed) where it is linked to the Basic UDI-DI. The URL to the Eudamed public website is: <https://ec.europa.eu/tools/eudamed>. The SSCP is also available online at remmed.com

CH

REP

MedEnvoy
Switzerland,
Gotthardstrasse 28,
6302 Zug,
Switzerland

EC

REP

EMERGO EUROPE
Westervoortsedijk 60
6827 AT Arnhem
The Netherlands

UK

REP

EMERGO Consulting (UK)
Limited c/o Cr360
- UL International
Compass House, Vision Park
Histon
Cambridge
CB24 9BZ

STERILE

EO

UK

CA

0086

CE

2797

MD

Rx Only

RoHS

COMPLIANT

LEAD-FREE

RoHS Compliant
Lead-Free

Manufacturer

Remington Medical, Inc.
6830 Meadowridge Court
Alpharetta, GA 30005 USA
+1.800.989.0057
quality@remmed.com

DO

2

STERILIZE

Do Not Re-sterilize

2

Do Not Re-use

Do Not Use if Packaging is Open or Damaged

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