
**301-CG**
DISPOSABLE EXTENSION CABLE

Basic UDI-DI: 081307902CABLEW7

INTENDED PURPOSE:

The Disposable Extension Cables are defined as surgically non-invasive devices intended for transient use. The Disposable Extension Cable Families are intended to be used by professional medical personnel during specific cardiac procedures such as pacemaker implantation and cardiac functionality analysis, etc. The cables are a single-patient use device and are not intended to be used for more than one patient nor for several procedures as these may increase risk of infection. The 301-CG Disposable Extension Cable is intended to act as a conduit for the electrical signal generated from a Pacing System Analyzer (PSA) to a permanent pacemaker lead. After connections are completed, analyzing, pacing, and sensing thresholds can be measured.

INDICATIONS FOR USE:

The 301-CG Disposable Extension Cables are electrical extension cables designed to transmit signal from patient-connected leads to a PSA. The 301-CG Disposable Extension Cables do not have any therapeutic or diagnostic indications or claims.

WARNINGS/PRECAUTIONS:

1. Transmittance of micro-shocks from static electricity can be generated through the pacing system. It is recommended that gloves be used when handling any part of the pacing system while in use.

2. Damage to cable may occur from improper handling.
3. DO NOT reuse. This may cause the cables to malfunction and/or to no longer fit with the PSA unit.
4. DO NOT use if sterile packaging is damaged or unintentionally opened before use.
5. DO NOT use for patients weighing less than 6.49 kilograms.

CONTRAINDICATIONS:

Use of Disposable Extension Cables is contraindicated where any pacing procedures are contraindicated (i.e., use in patients with cardiac conditions that contraindicate temporary or permanent pacing).

INSTRUCTIONS FOR USE:

Remove the 301-CG Disposable Extension Cable from the sterile package and connect it to the PSA or adaptor as follows:

1. Single-Chamber Pacing:

Requires one 301-CG Disposable Extension Cable. Connect the extension cable for the ventricular lead system to the ventricular connector block. Attach the alligator clips in either a unipolar or bipolar configuration.

a. *Single-Chamber Unipolar Leads*

Attach the black (negative) clip to the connector pin on the end of lead. Attach the red (positive) clip to the temporary indifferent electrode.

b. *Single-Chamber Bipolar Leads*

Attach the black (negative) clip to the connector pin on the end of lead. Attach the red (positive) clip to the positive pin or positive ring on the end of lead.

2. Dual-Chamber Pacing:

a. *Dual-Chamber Unipolar Leads*

Attach the black (negative) clip of the atrial cable to the connector pin on the end of the atrial lead and attach the black (negative) clip of the ventricular cable to the connector pin on the end of the ventricular lead. Attach the red (positive) clip of both cables to the temporary indifferent electrode.


b. *Dual-Chamber Bipolar Leads*

Attach the black (negative) clip of the atrial cable to the connector pin on the end of the atrial lead and attach the black (negative) clip of the ventricular cable to the connector pin on the end of the ventricular lead. Attach the red (positive) clip of atrial cable to the positive pin or positive ring on the end of the atrial lead and attach the red (positive) clip of the ventricular cable to the positive pin or positive ring on the end of the ventricular lead.

DISPOSAL:

1. Dispose of the 301-CG Disposable Extension Cable in an environmentally friendly way as contaminated medical waste.

Any serious incident which has occurred in relation to the 301-CG Disposable Extension Cable should be reported to Remington Medical, Inc. and to the competent authority of the Member State in which the user and/or patient is established. The Summary of Safety and Clinical Performance (SSCP) for Remington Medical Disposable Extension Cables is available in the European database on medical devices (Eudamed) where it is linked to the Basic UDI-DI. The URL to the Eudamed public website is: <https://ec.europa.eu/tools/eudamed>. The SSCP is also available online at remmed.com

Rx Only

Remington Medical, Inc.
6830 Meadowridge Court
Alpharetta, GA 30005 USA
+1.800.989.0057
quality@remmed.com

Do Not Re-sterilize

Do Not Re-use

Do Not Use if Packaging is Open or Damaged

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