DE

Remington MEDICAL

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Remington Medical, Inc.
6830 Meadowridge Court
Alpharetta, GA 30005 EUA
+1.800.989.0057
quality@remmed.com

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EN

Remington MEDICAL

AUTOMATIC CUTTING NEEDLE NAC Needle Family (NAC, NAC B, NAC M, NAC MLL, NAC U, NAC ULTRA) SOFT TISSUE BIOPSY NEEDLE

Indications for Use:

Automatic Cutting Needles are used to obtain multiple needle biopsies/core samples from soft tissue. These needles are to be used in conjunction with the following biopsy instruments as represented in the table below.

Biopsy Instrument:	Fits with Needle:
Angiotech® Pro-Mag™ 2.2	NAC, NAC B, NAC U
Angiotech® Pro-Mag™ ULTRA	NAC ULTRA
Bard® Biopsy®	NAC, NAC B, NAC U
Bard® MAGNUM®	NAC M, NAC MLL
MBD™-23	NAC, NAC B

Warnings/Precautions:

- 1 For use only for core biopsies of soft tissue.
- 2 Needle should be used by a physician familiar with the possible side effects, typical findings, limitations, indications and contraindications of a core needle biopsy procedure.
- 3 Physician judgment is required when considering a biopsy on patients with bleeding disorders, or receiving anticoagulant medications.
- 4 For single patient use only; if re-used, needle set may not correctly fit with applicable biopsy instrument or may bend/break during use.

Instructions for Use:

- 1 Open the sterile package using proper aseptic technique.
- 2 Remove the needle from the sterile package and remove the tip protector.
- 3 Inspect the needle (stylet and cannula) for damage that would prevent proper operation of the needle assembly
- 4 Open the lid of the biopsy device and load the needle set into the device.
- 5 Operate the biopsy device following the manufacture's operating instructions.

Note: After use, this product may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable laws and regulations. Handle in a manner which will prevent accidental needle puncture.

Sterilized by Ethylene Oxide

Do Not Reuse

Do Not Resterilize

Do Not Use if Package is Open or Damaged

Manufacture

Remington Medical, Inc.
6830 Meadowridge Court
Alpharetta, GA 30005 USA
+1.800.989.0057
quality@remmed.com

2797

EMERGO EUROPE
Prinsessegracht 20
2514 AP The Hague
The Netherlands

Rx Only

Non-Pyrogenic

Angiotech® and Pro-Mag™ are registered trademarks of Angiotech or its affiliates. Bard®, Biopsy® and MAGNUM® are registered trademarks of C.R. Bard, Inc. or its affiliates. MBD™ is a registered trademark of Promex Technologies, LLC or its affiliates.

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