


Remington
MEDICAL

PRODUCT CATALOG



Remington MEDICAL

At Remington Medical, we have a dedication to quality and control. That is why we have done our best to bring many manufacturing functions in-house. We believe that when we have control over all facets of the manufacturing process a superior product is the result. That is why we manufacture our raw components, assemble in-house, and pass the savings to you. Our goal is to reduce the cost to the end user, improve patient care, and ensure reliability.

Products include: pacing cables, needles, syringes, drainage bags and introducer kits.

At Remington Medical our goal is to provide safe disposable medical devices of the highest quality to our customers quickly with the greatest standard of customer service.

Founded in 1992 by Fred L. Aycock, Remington Medical has grown from a small 2,000 square feet office to an international medical device manufacturer. Remington Medical has grown to three facilities between Alpharetta, Georgia and Santo Domingo, Dominican Republic totaling more than 80,000 square feet.

With each great milestone, Remington Medical has continued to provide customers with high quality, cost effective products they can depend on consistently. Remington Medical manufactures and distributes disposable medical devices to several medical specialty fields. These include products for Cardiology, Cardiovascular Surgery, Oncology, Vascular Access, Urology and General Surgery.

By integrating all aspects of the product development cycle under one umbrella, Remington Medical is able to closely monitor every step taken in the design, development and completion of customer's products.

PRODUCT: CABLES

S-101-97



Disposable small and large alligator clip style cable with safe connect used for pacing procedures.

ADAP-2000



Disposable pacing cable with pac-loc safety device used for transvenous and temporary pacing.

FL-601-97




Disposable screw-down extension cable with remington safety separator for temporary pacing of open hearts.

BS-101-97



Disposable small alligator clip style cable with safe-connect used for bi-ventricular pacing procedures.



CABLES


800.989.0057 www.remmed.com

CABLES

800.989.0057 www.remmed.com

S-101-97

Disposable small and large alligator clip style cable with safe connect used for pacing procedures.




- Box of 5
- United States & Canada
- Sterile
- Large & Standard Sizes
- 8 ft, 12 ft
- Small and large alligator clips available for individual preferences
- Remington Safe-Connect for a secure fit to the EPG or PSA
- Compatible with Medtronic 5348 or 5388 EPG
- Safe, cost effective and reliable alternative to reusable pacing cables



- Box of 5
- United States & Canada
- Sterile
- The Pace-Loc Clip Secures the connection between the catheter and adapter
- Compliance with FDA standards
Patient safety is not compromised
- Meets industry IEC 60601-1 standard for design, functionality and safety
- 6 ft, 12 ft

ADAP-2000

Disposable pacing cable with Pac-Loc safety device used to transvenous and temporary pacing. The ADAP-2000 cable is specifically designed to provide a safe connection between 2mm shrouded pin pacing catheters and the Medtronic Temporary Pacemaker Model 5348 or 5388 (or Medtronic 5101 Adapter and ADAP - I Adapter).



CABLES


800.989.0057 www.remmed.com

CABLES

800.989.0057 www.remmed.com

FL-601-97

Disposable screw-down extension cable with Remington Safety Separator for temporary pacing of open hearts. The FL-601-97 Fas-Loc accepts an epicardial wire at one end and plugs into a Medtronic 5348 or 5388 EPG or older EPG with an adapter (Medtronic 5101 or Remington ADAP I Adapter).



- Box of 5
- United States & Canada
- Sterile
- Atrial, Atrial/Ventricular
- 6 ft, 12 ft

Remington safety separator (U.S. Patent No. 7059878 & 7204703) designed to ensure the epicardial wires remain separate

Remington safe-connect for a secure fit to the EPG

Safe, cost effective, and reliable alternative to reusable pacing cables



- Box of 5
- United States & Canada
- Sterile

Small and large alligator clips available for individual preferences


Remington Safe-Connect for a secure fit to the EPG or PSA

Compatible with Medtronic 5348 or 5388 EPG

Safe, cost effective and reliable alternative to reusable pacing cables

BS-101-97

Disposable small alligator clip style cable with safe-connect used for bi-ventricular pacing procedures.



PRODUCT: NEEDLES & GUIDES

Bone Marrow Aspiration Needle



Bone Marrow Aspiration Needle with Adjustable Needle Stop. Safe and Simple Aspiration from the Sternum or Iliac Crest.

Bone Marrow Biopsy Needle



Unique "Double Diamond" Stylet Point is Designed for Quick and Easy Penetration into the Marrow Cavity. Available in a Variety of Sizes

Breast Localization Needle



Needle Designed to Pinpoint a Suspicious Lesion Within the Breast Tissue to Be Biopsied (Note: Wire Is Intended for One-Time Placement only - Surgical Removal Required).

All Style Needles



Disposable biopsy needles used to obtain core tissue samples for cytological evaluation. Compatible with the Bard® Magnum® biopsy instrument.



NEEDLES

800.989.0057 www.remmed.com

NEEDLES

800.989.0057 www.remmed.com

Bone Marrow Aspiration Needle

Bone Marrow Aspiration Needle With Adjustable Needle Stop. Safe And Simple Aspiration From The Sternum Or Iliac Crest.



Box of 10

United States & Canada

Sterile

15 Gauge, 16 Gauge

1 7/8 IN

Ergonomically designed handle combined with the match ground cannula and stylet points provide clinician comfort while easily penetrating hard bone

Safe aspiration at the sternum can be performed by adjusting the needle stop to control the depth of penetration

Needle stop can be removed if aspirating at the iliac crest

Luer lock connector on the handle provides for secure syringe attachment

Available in two-gauge sizes



Box of 10

United States & Canada

Sterile

11 Gauge

4 IN, 6 IN

Ergonomically designed handle allows clinician to apply sufficient pressure for penetration without excessive hand fatigue


Luer lock connector on the handle provides for secure syringe attachment

Needle stop can be removed if aspirating at the iliac crest

Available in a variety of sizes

Bone Marrow Biopsy Needle

Unique "Double Diamond" stylet point is designed for quick and easy penetration into the marrow cavity.



NEEDLES

800.989.0057 www.remmed.com

NEEDLES

800.989.0057 www.remmed.com

Breast Localization-Needle

Needle designed to pinpoint a suspicious lesion within the breast tissue to be biopsied. (Note: wire is intended for one-time placement only - surgical removal required).



- Box of 10
- United States & Canada
- Sterile
- 20 Gauge
- 3 CM, 5 CM, 7 CM, 10 CM
- Includes one wire, which is pre-loaded for convenience.
- Kopans-style needle
- Easy thread wire entry
- Stiffened wire pre-loaded in needle



- Box of 10
- United States & Canada
- Sterile
- 18 Gauge
- 16 CM, 20 CM, 25 CM
- Universal design fits most biopsy devices including Bard® Biopty® & Argon Pro-Mag™ 2.2
- Provides large & fully intact specimens
- Echogenic enhanced tip for placement under ultrasound guidance
- Centimeter markings on the outer cannula to facilitate depth placement
- Color coded hubs indicate gauge size

B-Needle Soft Tissue Biopsy Needle

Disposable biopsy needle used to obtain core tissue samples for cytological evaluation. Fits most biopsy instruments including Bard® Biopty® & Argon Pro-Mag™ 2.2.

NEEDLES

800.989.0057 www.remmed.com



Chiba Style Needle

General purpose aspiration and injection needle.



Box of 10

United States & Canada

Sterile


18,20,22,23 Gauge

9 CM, 15 CM, 18 CM, 20 CM, 25 CM

- Short Bevel Removes Sample From Controlled Area
- Centimeter Markings On The Outer Cannula To Facilitate Precise Depth Placement
- Adjustable Needle Stop Restricts Forward Movement Localizing The Needle Tip To The Biopsy Site

NEEDLES

800.989.0057 www.remmed.com



BB-Needle Soft Tissue Biopsy



Box of 10


United States & Canada

Sterile

18 Gauge

16 CM, 20 CM, 25 CM

- Recommended For Use With The Bard® Biopty® Biopsy Instrument or Argon Pro-Mag 2.2.
- Provides Large & Fully Intact Specimens
- Echogenic Enhanced Tip For Placement Under Ultrasound Guidance
- Centimeter Markings On The Outer Cannula To Facilitate Depth Placement
- Color-Coded Hubs Indicate Gauge Size



NEEDLES

800.989.0057 www.remmed.com

NEEDLES

800.989.0057 www.remmed.com

M-Needle Soft Tissue Biopsy Needle

Disposable biopsy needle used to obtain core tissue samples for cytological evaluation. Compatible with the Bard® Magnum®.



- Box of 10
- United States & Canada
- Sterile
- 18 Gauge
- 16 CM, 20 CM, 25 CM
- Recommended for use with the Bard® Magnum®
- Provides large & fully intact specimens
- Echogenic enhanced tip for placement under ultrasound guidance
- Centimeter markings on the outer cannula to facilitate depth placement
- Color coded hubs indicate gauge size



- Box of 10
- United States & Canada
- Sterile
- 14, 18 Gauge
- 16 CM, 20 CM, 25 CM
- Recommended for use with the Bard® Magnum® biopsy instrument
- Provides large & fully intact specimens
- Echogenic enhanced tip for placement under ultrasound guidance
- Centimeter markings on the outer cannula to facilitate depth placement
- Cannula hub accepts a luer lock syringe for aspirations & injections
- Color coded hubs indicate gauge size

MLL-Needle Soft Tissue Biopsy Needle

Disposable biopsy needle used to obtain core tissue samples for cytological evaluation. Compatible with the Bard® Magnum® biopsy instrument.

NEEDLES

800.989.0057

www.remmed.com

NEEDLES

800.989.0057

www.remmed.com

N-Style Soft Tissue Biopsy Needle

Disposable Biopsy Needle Used To Obtain Core Tissue Samples For Cytological Evaluation. Compatible With The Argon Pro-Mag™ 2.2 Biopsy Instrument.



Franseen Style Biopsy Needle

General Purpose Aspiration and Injection Needle.

Box of 10

United States & Canada

Sterile

14, 18 Gauge

10 CM, 16 CM, 20 CM

- Recommended For Use With The Argon Pro-Mag™ 2.2 Biopsy Device
- Echogenic Enhanced Tip For Placement Under Ultrasound Guidance
- Centimeter Markings On The Outer Cannula To Facilitate Depth Placement
- Color-Coded Hubs Indicate Gauge Size

Box of 10


United States & Canada

Sterile

18 , 20 , 22 Gauge

18 CM, 20 CM, 22 CM

- Cannula Point Geometry Suitable For Cutting Difficult Lesions
- Centimeter Markings On The Outer Cannula To Facilitate Precise Depth Placement
- Adjustable Needle Stop Restricts Forward Movement Localizing The Needle Tip To The Biopsy Site



NEEDLES

800.989.0057 www.remmed.com

NEEDLES

800.989.0057 www.remmed.com

U-Style Soft Tissue Biopsy Needle

Disposable Biopsy Needle Used To Obtain Core Tissue Samples For Cytological Evaluation. Universal Design Fits Most Biopsy Devices Including Bard® Biopty® & Argon Pro-Mag™ 2.2.



Ultra-Style Soft Tissue Biopsy Needle

Disposable Biopsy Needle Used To Obtain Core Tissue Samples For Cytological Evaluation. Compatible With The Argon Pro-Mag Ultra 2.5 Biopsy Instrument.

Box of 10

United States & Canada

Sterile

18 Gauge

20 CM, 25 CM

Universal Design Fits Most Biopsy Devices Including Bard® Biopty® & Argon Pro-Mag™ 2.2

Provides Large & Fully Intact Specimens

Echogenic Enhanced Tip For Placement Under Ultrasound Guidance

Centimeter Markings On The Outer Cannula To Facilitate Depth Placement

Color-Coded Hubs Indicate Gauge Size

Box of 10

United States & Canada

Sterile

18 Gauge

20 CM, 25 CM

Recommended For Use With The Argon Pro-Mag Ultra 2.5 Biopsy Device

Echogenic Enhanced Tip For Placement Under Ultrasound Guidance

Centimeter Markings On The Outer Cannula To Facilitate Depth Placement

Color Coded Hubs Indicate Gauge Size

REMININGTON

PRODUCT: SUPPLY

600 ml
Drainage
Bag



Disposable 600ml Drainage Bag Ideal For Fluid Drainage Including Abscess, Biliary, And Nephrostomy.

Locking
Syringe



10 cc Locking Syringe

Bag
Decanter




Spiked Bag Decanter for the Transfer of Sterile Fluids.

Tear Away
Sheath
Introducer
Kit



Introducer Set Used For Placement Of A .038 in. Diameter Guide Wire Into The Vascular System



600 ml Drainage Bag

Disposable 600ml Drainage Bag Ideal For Fluid Drainage Including Abscess, Biliary, And Nephrostomy.



Box of 20

United States & Canada

Sterile

- 24 Soft, Kink-Resistant Drainage Tube Can Be Cut To Appropriate Length
- Anti-Reflux Valve To Prevent Retrograde Contamination
- Fully Adjustable Soft Waist Belt
- Twist Type Outlet Valve For Quick And Easy Emptying
- Not Made Of Natural Rubber Latex
- Single Use



Box of 10

United States & Canada

Sterile


Angle: 15 Degree, 35 Degree

- Not Made Of Natural Rubber Latex
- Choices Of Angles For Ease Of Dispensing
- Includes Two Caps
- One Piece, Rigid Construction
- Transparent For Visual Flow
- Packaged In Sterile Pouch
- The 9001 Bag Decanter Features a 15 Angle
- The 9002 Bag Decanter Features a 35 Angle

Bag Decanter

Spiked Bag Decanter for the Transfer of Sterile Fluids.





Locking Syringe

10 cc Locking Syringe



Box of 25

United States & Canada

Sterile

10 cc Luer Lock Syringe With Locking Mechanism
Locks In At 6 cc And Holds Vacuum For Aspiration Biopsy Procedures

Box of 10

United States

Sterile

Introducer Size:
6,7,8,9,10,11 French



Tear Away Sheath Introducer Kit

Introducer Set Used For Placement Of A .038 in. Diameter Guide Wire Into The Vascular System.

Kink Resistant PTFE Sheath Ensures an Even Tear Every Time

Tear away introducer sheath with Spin-Lock For Ease Of Insertion

Smooth Dilator Taper And Dilator-To-Sheath Transition

Custom Sizes Available

Complete Introducer Kit Includes:

Dilator/Tear-Away Sheath Assembly

18 Gauge Arterial Needle

12 cc Syringe

.038 in. x 50 cm Guide Wire

Automatic Cutting Biopsy Needles (NAC)

TYPE	Available Lengths				Units Available			Gauge Size
	10 CM	16CM	20CM	25CM	Sterile		Non-Sterile	
					Box (10 Units)	Case (10 boxes)	Bag	
N-Needle	NAC 1410				✓	✓		14 GA
		NAC-1816	NAC-1820		✓	✓		18 GA
B-Needle					✓	✓		14 GA
		NAC-1816B	NAC-1820B	NAC-1825B	✓	✓		18 GA
BB-Needle					✓	✓		14 GA
		NAC-1816BB	NAC-1820BB	NAC-1825BB	✓	✓		18 GA
M-Needle					✓	✓		14 GA
		NAC-1816M	NAC-1820M	NAC-1825M	✓	✓		18 GA
MLL-Needle	NAC-1410MLL				✓	✓		14 GA
			NAC-1820MLL	NAC-1825MLL	✓	✓		18 GA
U-Needle					✓	✓		14 GA
			NAC-1820U	NAC-1825U	✓	✓		18 GA
Ultra Needle					✓	✓		14 GA
			NAC-1820 ULTRA	NAC-1825 ULTRA	✓	✓		18 GA

Bone Marrow Needles

TYPE	Available Lengths			Units Available			Gauge Size
	1 7/8 IN	4 IN	6 IN	Sterile		Non-Sterile	
				Box (10 Units)	Case (10 boxes)	Bag (25 Units)	
Aspiration Needle	1501-I			✓	✓		15 GA
	1601-I			✓	✓		16 GA
Double Diamond Needle		JD-1104	JD-1106	✓	✓	✓	11 GA

Breast Localization Needles

TYPE	Available Lengths				Units Available			Gauge Size
	3 CM	5 CM	7 CM	10 CM	Sterile		Non-Sterile	
					Box (10 Units)	Case (10 boxes)	Bag	
Breast Localization Needles	BNL-2003	BNL-2005	BNL-2007	BNL-2010	✓	✓		20 GA

Chiba Needles

TYPE	Available Lengths					Units Available			Gauge Size
	9 CM	15 CM	18 CM	20 CM	25 CM	Sterile		Non-Sterile	
						Box (10 Units)	Case (10 boxes)	Bag (125 Units)	
Chiba Needle	CNM-1803	CNM-1806		CNM-1808		✓	✓	✓	18 GA
	CNM-2003	CNM-2006		CNM-2008		✓	✓	✓	20 GA
	CNM-2203	CNM-2206	CNM-2007	CNM-2208	CNM-2210	✓	✓	✓	22 GA
		CNM-2306		CNM-2308		✓	✓	✓	23 GA

Franseen Needles

TYPE	Available Lengths			Units Available			Gauge Size
	9 CM	15 CM	18 CM	Sterile		Non-Sterile	
				Box (10 Units)	Case (10 boxes)	Bag (125 Units)	
Franseen Needle	FNM-1806	CNM-1808		✓	✓		18 GA
	FNM-2003	FNM-2006	FNM-2008	✓	✓		20 GA
	FNM-2203	FNM-2206	FNM-2008	✓	✓		22 GA

Guidewire Introducer Needles

TYPE	Available Lengths	Units Available			Gauge Size
	7 CM	Sterile		Non-Sterile	
		Box (10 Units)	Case (10 boxes)	Bag (200 Units)	
Standard Round Hub	GW1-1802	✓	✓	✓	18 GA
	GW1-1902	✓	✓	✓	19 GA
Standard Square Hub	GW1-1802WW	✓	✓	✓	18 GA
	GW1-1902WW	✓	✓	✓	19 GA
Funnel Tapered Hub with Wings	BGWI-1802	✓	✓	✓	18 GA
	BGWI-1902	✓	✓	✓	19 GA

TYPE	Available Lengths				Units Available			Gauge Size
	2 CM	3 CM	5 CM	6 CM	Sterile		Non-Sterile	
					Box (10 Units)	Case (10 boxes)	Bag (125 Units)	
N-Needle		TNMR-1703	TNMR-1705		✓	✓	✓	17 GA
		TNMR-1803	TNRM-1805		✓	✓	✓	18 GA
	TNMR-1902				✓	✓	✓	19 GA
		TNMR-2003	TNMR-2005	TNMR-2006	✓	✓	✓	20 GA
		TNRM-2203			✓	✓	✓	22 GA
M-Needle		TNMF-1703	TNMF-1705		✓	✓	✓	17 GA
		TNMF-1803	TNMF-1805		✓	✓	✓	18 GA
					✓	✓	✓	19 GA
	TNMF-1902	TNMF-2003	TNMF-2005	TNM-2006	✓	✓	✓	20 GA
		TNMF-2203			✓	✓	✓	22 GA

Westcott Needles

TYPE	Available Lengths		Units Available			Gauge Size
	7 CM	15 CM	Sterile		Non-Sterile	
			Box (10 Units)	Case (10 boxes)	Bag	
Westcott Needle	WNM-2203	WNM-2006	✓	✓		20 GA
	WNM-2203	WNM-2006	✓	✓		22 GA

Cable Adapters

Product Code	Description	Units	Attribute
ADAP-2R	Redel Connector Adapter	1 ea	Non-sterile
ADAP-I	Cable Adapter	1 ea	Non-sterile
ADAP-II	Lemo Connector Adapter	1 ea	Non-sterile

TYPE	Available Lengths			Units Available			Variant
	6 FT	8 FT	12 FT	Sterile		Non-Sterile	
				Box (5 Units)	Case (20 boxes)	Bag (25 Units)	
Open Heart Cable	FL-601-97		FL-601-97-12	✓	✓	✓	Atrial
	FL-601-97-B		FL-601-97-12-B	✓	✓	✓	Atrial/Ventricular
Small Alligator Clip		S-101-97	S-101-97-12	✓	✓	✓	Atrial
		S-101-97-B		✓	✓	✓	Atrial/Ventricular
Large Alligator Clip		S-201-97	S-201-97-12	✓	✓	✓	Atrial
		S-201-97-B		✓	✓	✓	Atrial/Ventricular
Biventricular Small Alligator Clip		BS-101-97	BS-101-97-12	✓	✓	✓	Atrial
			NAC-1820M	✓	✓	✓	Atrial/Ventricular
Shrouded Pins		PA-801		✓	✓	✓	Atrial
				✓	✓	✓	Atrial/Ventricular
Biotronik Connector				✓	✓	✓	Atrial
							Atrial/Ventricular
Temporary Cable							
	ADAP-2000		ADAP-2000-12				

Introducer Kits

Product Code	Size	Units Available Sterile Box (10 units)
7006	6 FR	✓
7007	7 FR	✓
7008	8 FR	✓
7009	9 FR	✓
70010	10 FR	✓
70011	11 FR	✓

Bag Decanters

		Units Available		
Product Code	Angle (degrees bent)	Sterile		Non-Sterile
		Box (50 Units)	Case	Bag (250 Units)
9001	35	✓		✓
9002	35	✓		✓

Drainage Bags

		Units Available		
Product Code	Size	Sterile		Non-Sterile
		Box (20 Units)	Case	
600-D	600 ml	✓		

Locking Syringes

		Units Available		
Product Code	Size	Sterile		Non-Sterile
		Box (25 Units)	Case	
LSS-100	10 cc	✓		

Needle Guides

Product Code	Attributes	Available Lengths			Units Available		
		14 GA	20 GA	22 GA	Sterile		Non-Sterile
					Box (10 Units)	Case	Bag
LNGC-14		✓	✓	✓	✓	✓	
					✓		
M-Needle	Closed		✓	✓	✓		
	Open						